



Faculty of Electrical and Computer Engineering - Institute of Semiconductors and Microsystems

FTIR-ATR spectroscopy with novel microstructured single reflection elements (mSRE) basing on silicon wafers as internal reflection elements

Ulrich Künzelmann

nanoSeminar

tud » faculty of mechanical engineering » institute for materials science » chair "materials science and nanotechnology" Dresden, August 25th, 2011



Introduction/Motivation

- Novel microstructured single reflection elements at silicon wafers for ATR-FTIR spectroscopy
- In situ-Investigations
 - Wet etch of thermal SiO₂ layers at Si
 - Chemical mechanical planarisation of SiO₂ layers
 - Native SiO₂ layers and surface termination at Si
 - Interactions of additives for Cu electrodeposition
 - General applications in (bio)analytical chemistry
- Conclusion / Outlook

Motivation



Motivation | Novel Technique | In situ Investigations | Conclusions

- Characterisation of thin-films on Si/SiO₂ substrates
- Original focus on chemical-mechanical polishing (CMP) and electrochemical deposition (ECD):
 - Wet chemical processes (e.g. with aqueous solutions)
 - Surface reactions
 - Tribochemical/physicochemical investigations

 \rightarrow In situ measurement (

Surface not accessible with common surface analytical methods like XPS, XRD, SIMS, LEED (mostly vacuum at interesting face)

- ► Fourier transform infrared (FTIR) spectroscopy
 - Molecular spectroscopy (information on composition and structure)
 - Based on change in dipole moment due to molecule vibrations
 - Vacuum not necessary, but applicable at each face, if advantageous

TECHNISCHE UNIVERSITAT Internal vs. external reflection technique



Motivation | Novel Technique | In situ Investigations | Conclusions

- External reflection: IR penetrates ambient and layer, IRRAS
- Internal reflection: IR penetrates substrate and is attanuated by layer, ATR
- Advantages of the ATR technique:
 - Measurement of thick or strong absorbing media
 - No interferences on thin samples/layers
 - Surface sensitivity
 - Works with polarisation



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 4/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Microstructured SRE for ATR-FTIR



Motivation | Novel Technique | In situ Investigations | Conclusions

- ► Common Si {100} wafer
- Microstructure on the wafer backside
 - V-shaped grooves of Si(111) faces
 - common MEMS fabrication technologies (e.g. wet chemical crystal oriented etch)
- IR beam from the wafers backside, optical path can be purged or evacuated
- Attenuated total reflection (ATR) at the Si/SiO₂ or Si/ambient interface
- Penetration depth can be partially controlled
- Horizontal or vertical sample orientation
- Sample side is free accessible



H. Schumacher, U. Künzelmann, B. Vasilev, J. W. Bartha, and K. J. Eichhorn "Applications of Microstructured Silicon Wafers as Internal Reflection Elements in ATR-FTIR Spectroscopy" *Applied Spectroscopy* **64** (9), 2010, p. 1022-1027

nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 5/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de



Motivation | Novel Technique | In situ Investigations | Conclusions

 mSRE has a very short optical path (≈1mm) in Si (70mm with MRE)
 No limitation due to Si lattice vibrations
 Enhanced measurement range for Si reflection elements (entire mid and far infrared spectral range)



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 6/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Wet etch of thick thermal oxide



Now accessible with mSRE!

Motivation | Novel Technique | In situ Investigations | Conclusions

- ▶ Wet etch of SiO₂ with buffered oxide etch solution (H₂O, HF, NH₄F)
- ▶ *In situ* ATR-FTIR (DTGS detector, 4 cm⁻¹, 64 scans)
- Decrease of SiO₂ assigned peaks
- Increase of BOE assigned peaks (e.g. H₂O)





U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 7/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Experimental setup (CMP)



Motivation | Novel Technique | In situ Investigations | Conclusions

- Single reflection measurement cell
- Angle of incidence of 35° on {100} wafer backside
- ▶ Information depth
 ≈ 900 nm @ 1000 cm⁻¹





Schematic of the polishing configuration and of the simulator at the FTIR instrument

nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 8/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de In situ measurement during CMP



Motivation | Novel Technique | In situ Investigations | Conclusions

- ► Silica based, particle size ≈50 nm; 30 wt.% solids; pH 10; NH₄OH chemistry
- Pad: IC 1000, conditioned, linear movement
- ► Back pressure ≈ 11 psi
- ▶ FTIR: DTGS detector, 64 scans/spectrum, 4 cm⁻¹ spectral resolution



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 9/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

In situ measurement during CMP



Motivation | Novel Technique | In situ Investigations | Conclusions



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 10/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Impact of the polishing pad



Motivation | Novel Technique | In situ Investigations | Conclusions

- Not observed during polishing
- Observable during static measurement



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 11/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Thin layers and surface termination



Motivation | Novel Technique | In situ Investigations | Conclusions

- ▶ Wet chemical etch of an approx. 2 nm thin SiO₂ with diluted HF-solution
- ► Native oxide is removed → *negative peak*
- When reaching bulk silicon Si-H vibrations arise \rightarrow positive peak

 \rightarrow Surface termination with hydrogen (silane groups)



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 12/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Interaction of plating additives



Motivation | Novel Technique | In situ Investigations | Conclusions

- Suppressor (e.g. *polyethylene glycol*, *PEG*) adsorbs at exposed areas and inhibits the Cu deposition → global polarisation
- ► Accelerator (e.g. *bis-(3-sulfopropyl)-disulfide, SPS*) displaces suppressor at the bottom (forming (3-mercaptopropyl)sulfonate) → local depolarisation



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011 U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 13/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de



Motivation | Novel Technique | In situ Investigations | Conclusions

In-situ-spectroelectrochemical measurements during the reduction of methylviologene



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 14/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Application in spectroelectrochemisty II

Motivation | Novel Technique | In situ Investigations | Conclusions



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 15/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Intrinsically conductive polymers



Motivation | Novel Technique | In situ Investigations | Conclusions



nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011 U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 16/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Intrinsically conductive polymers



Motivation | Novel Technique | In situ Investigations | Conclusions



nanoSeminar, Chair of "Materials Science and Nanotechnology" U. K TUD/Faculty Mechanical Engineerung/Inst. of Material Science ulric Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 17/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de





nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 18/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de

Conclusion



Motivation | Novel Technique | In situ Investigations | Conclusions

- Microstructured SRE for ased on backside structured Si-wafers
 - \rightarrow No limitations in sample size
- ► ATR-FTIR analysis from the wafer backside
- Enhanced usable spectral range for Si based elements
 - \rightarrow No limitations due silicon lattice vibrations
 - \rightarrow below 1500 cm⁻¹
- ▶ In situ investigations of wet chemical and *tribochemical* processes
- Detection and analysis of
 - slurry incl. abrasive and additives
 - (ultra) thin-films
 - surface allocations
 - additives for Cu electroplating
- Standard applications for FTIR measurements with cheap RE





- Microstructured SRE applicable for investigation of polymers
- ► Biological samples can be measured → As bulk substances
 - \rightarrow In the encapsulated/immobilized state
- Application of barrier and/or diffusion layers provides surface selectivity
- Antibody-antigene interactions observable
- Micro-integration possible

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 20/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de



Best thanks to:

- H. Schumacher and Prof. J. W. Bartha, IHM.
- The processing team of the *IHM* clean room.
- G. Adam and Dr. K.-J. Eichhorn from the *Leibniz-Institut für Polymerforschung – IPF Dresden* for proficiently supporting the infrared measurements in the initial phase of the work.
- Dr. K. Hinrichs from the *Leibniz–Institut für Analytische Wissenschaften - ISAS Berlin* for helpful discussions.
- The financial support by the *Bundesministerium für Bildung und Forschung* (BMBF, FKZ: 01M3183B and 13N10808).

U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 21/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de





nanoSeminar, Chair of "Materials Science and Nanotechnology" TUD/Faculty Mechanical Engineerung/Inst. of Material Science Dresden, August 25th, 2011 U. Künzelmann.: FTIR-ATR spectroscopy with novel mSRE ... 22/20 ulrich.kuenzelmann@tu-dresden.de www.ihm.tu-dresden.de